FR

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Signature & Company stamp

EN

Bluetooth Body Fat Scale
Mode: FG210
INSTRUCTION MANUAL

Read this manual thoroughly before using and save it for future reference

FitGo - Ideas for health

Health & fitness

iOS App

Android App

One-click sync with

Mobile&Device

iOS 7.0

Android 4.3

Bluetooth 4.0

For First Use

1. Battery Installation

Open the battery door, put into batteries with correct polarity.

2. App Installation

Scan the QR code on the above, or download "FitGo" from App Store or Android Google Play.

3. Registration

Turn on bluetooth on mobile device, run the App "FitGo", follow the App instructions.

12 parameters in 1 measurement

Weight

BMI

Body Fat

Water

Muscle

Bone

BMR

Obesity

V-Fat

Protein

AGE

LBM

Multi Users

8 users share one user account in the device.

Target weight

Users can set their own target weight to customize their own health plan.

Attention

1. Do not use the scale on people who have body implants such as a pacemaker, artificial limbs, contraceptive devices, or metal plates or screws. It may cause the device to malfunction or produce an inaccurate result. When in doubt, consult your physician.

2. Please make sure good contact between bare feet and the metal electrodes, otherwise it will display your weight only.

3. Ensure the scale is on a flat, even and firm surface. Avoid carpet or soft surfaces.

4. The surface of the scale will be slippery if wet. Please make sure both the glass and your feet are dry before use.

5. This scale is not water proof. Do not expose to moisture.

6. Clean the scale with a slightly damp cloth. DO not use any chemical/abrasive cleaning agents to clean the scale.

7. Please remove the batteries, if the scale is not to being

used for a long time.

8.The appliance is intended solely for domestic use within the home.

9. Children should be supervised to ensure that they do not play with the appliance.

10. Frequency range: 2402MHz ~ 2480MHz

Battery Attention:

- Contains 2 AAA batteries, which can cause serious burns if swallowed, potentially resulting in internal bleeding or death. Keep out of sight and out of reach of small children and ensure battery compartment is screwed shut.
- Whenever you change the batteries make sure your hands are dry.
- Insert the batteries AAA in the battery compartment, respecting the polarity. The battery can only be replaced by batteries of the same type or equivalent. Replace all batteries at once. Do not mix new and old batteries.
- Remove all the batteries when the appliance will not be used for a long period of time and store them in a cool and dry place.
- Non-rechargeable batteries are not to be recharged.
- Different types of batteries or new and used batteries are not to be mixed.
- Batteries are to be inserted with the correct polarity.
- Exhausted batteries are to be removed from the appliance and safely disposed of.
- If the appliance is to be stored unused for a long period, the batteries should be removed.
- The supply terminals are not to be short-circuited.

CAUTION :

Risk of fire or explosion if the battery is replaced by an incorrect type.
Replace only with the same type.

REMOVAL OF APPLIANCES USED

Directive 2012/19/EU on waste electrical and electronic Equipment (WEEE) requires that the appliance uses are not disposed of using the normal municipal waste. Old appliances must be collected separately in order to optimize the recovery and recycling of materials they contain and reduce the impact on human health and the environment.

The symbol of the dustbin is used on all products to remind the obligations for separate collection. Consumers should contact their local authorities or dealer regarding the steps to follow for the removal of their old appliance. If you proceed with the scrapping of old equipment, be sure to render useless what could be dangerous: disconnect the power cable flush with the device.

Danger:

Never attempt to use your appliance, if it shows signs of damage or the power cord or plug is damaged. If the power cord is defective it is important to have it replaced with the after sales service to eliminate any danger.

Only qualified and authorized persons are authorized to carry out repairs to your appliance. Any repairs not in accordance with standards could significantly increase the level of risk for the user!

Defect resulting from improper handling, degradation or attempted repairs by third parties voids the warranty on the product. This also applies in case of normal wear and accessories from the appliance.

Important!

We recommend that you keep the packaging of your device at least for the duration of the warranty. The guarantee does not apply without the original packaging.

Warranty:

Your device has two years warranty. The warranty does not cover wear and breakage following a wrong use of the product. The customer is responsible for all returns. You remain responsible for the costs and risks associated with the product shipment; it is therefore recommended to ship the product with an acknowledgement of receipt and transport insurance in case of significant value of the product.

KEEP OUT OF REACH OF CHILDREN

Children should be supervised to ensure that they do not play with the appliance.

Service E-mail : relations.clients@euromenage.com

Declaration of Conformity

We, Zhongshan Canny Electronics Co., Ltd., declare under our sole responsibility that the product: FG210 is in conformity with the essential requirements and other relevant requirements of the RED Directive (2014/53/EU)

The complete CE declaration is available at the following address: <http://www.sgs.com/en/Terms-andConditions.aspx>

EU DECLARATION OF CONFORMITY

Manufacturer / Importer certifies that the following designated product

Product Description: Bluetooth Scale
Model No.: CF370BT (EUROMENAGE's model no. FG210/FG629)
Trade Mark: FAGOR

Is in conformity with the standard(s)

EN 55014-1:2017+A11:2020 EN55014-2:2015 EN301489-1
V2.2.3; EN301489-17 V3.1.1 EN300328 V2.2.2 EN62479:2010

Under the Council Directive

2014/53/EU on Electromagnetic Compatibility

The designated product was subject to sample testing for which
SGS issued test report number: GZEM210100003401,
GZEM210100003402, GZEM210100003403, GZEM210100003404

The declaration is the sole responsibility of the manufacturer / importer

Zhongshan Canny Electronics Co., Ltd.

2021-3-29

Date

Signature & Company stamp

EUROMENAGE
20/22 Rue de la ferme saint ladre
95470 SAINT WITZ FRANCE